RT-1300

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RT-1300

Manual of atmosphere annealing furnace

Contents

一、 Overview	3
二、 Technical parameter	3
三、 Furnace structure	3
四、 Intelligent temperature controller 518P	4
1.Introduction of intelligent temperature controller	4
2.Display switching of intelligent temperature controller	4
3.Program setting of intelligent temperature controller	6
4.Running of intelligent Temperature controller	7
5.Pause of the intelligent temperature controller	1
6.Stop of the intelligent temperature controller	1
五、 Maintenance Precautions	1

Thank you for purchasing the RT-1300 atmosphere annealing furnace. In order to prevent damage to the experimental furnace due to misuse, please read this instruction manual carefully before use.


一、 Overview

The RT-1300 atmosphere annealing furnace uses high-quality silicon carbide rods as heating elements, and the maximum furnace temperature can reach 1300°C. 30-segment heating curve can be programmed to heat up, automatic operation, and good temperature tracking. The electric furnace is widely used in the research of metal materials, ceramics, glass, etc.

二、 Technical parameter

name	unit	parameter
power	KW	9
voltage	V HZ	220V, 50HZ
maximum temperature	℃	1300
normal working temperature	℃	1200
heating rate	℃/min	<20
control precision	℃	±1℃
Heating element	Heating wire	Silicon carbide rod
Thermocouple graduation number	graduation	S type
Furnace size	mm	300*200*200mm

三、 Furnace structure

- | | | | |
|----------------------------------|------------------------------|------------------|------------------------------|
| 1 temperature control instrument | 5 power indicator light | 9 flow meter | 13 exhaust switch |
| 2 ammeter | 6 start button | 10 intake switch | 14 pressure gauge |
| 3 Voltmeter | 7 Stop button | 11 Exhaust port | 15 Pressure regulating valve |
| 4 Fault indicator light | 8 Control panel power switch | 12 Air inlet | 16 Vacuum pump socket |

4.0、 Introduction of intelligent temperature controller 518P

4.1.Main features:

(1) Adopt advanced AI artificial intelligence adjustment algorithm, no overshoot, with self-tuning function, can realize temperature rise and fall control with any slope, and have operation commands such as jump (cycle), run, pause and stop. Measurement accuracy: 0.25 level.

(2) 30 segments PID

temperature control


program control

function.

(3) Save data when

power off

2、 Panels description

(1) Furnace Temperature (P V)

(2) Target Temperature (S V)

(3) Set key (Confirm key) 

(4) Data shift key(Program setup entry)  (A/M)

(5) Data down key(Run/Hold too)  (RUN/HOLD)

(6) Data up key(Stop too)  (STOP)

(7) Function indication


(8) Percentage of power output (Light half stands for 50% output of rating power, lightning full standing for 100% output of rate power)

2.Display switch of intelligent controller

The working display of digital display indicates the working state of the controller, and its working state determines whether you can perform certain operations, so the user should pay attention to the working status of the controller panel when using the device or performing certain operations.

1、 Boot state:

After a few seconds of

temperature controller model
No. and software version
displayed, the basic state of
temperature measurement
will be displayed. "SV"
flashing and displays
"STOP", indicating that
the program is in a stopped
state, as shown in the figure
at right side

Model No. and Version


Basic state

2、Display switch as below instruction:

(1) In the basic state or program running state, press the XX key for 1 second to switch to the interface display (PV-STEP, SV-XX segment) programs running segment state. (Set the running segment or display the running temperature segment)

(2) Press XX the key again for 1 second to switch to the state of the running time of the segment. (Display the total running time of the running segment PV xxxx minutes, and elapsed running time SV xxxx minutes correspondingly)

(3) Press XX the key again for 1 second to return to the basic state.

3. Programs setting of intelligent temperature controller

The setting of the temperature control program is the user's choice of the process conditions of the sintering material, and the correct setting of the temperature control program is the prerequisite for the successful sintering of the material.


1. Programs setting as below:

(1) under basic state, Press YY for 1 second, controller enter the programs setting state. Controller firstly displayed the given value of the the present running segments, modify the value by press XYZ. The beginning value is 0 generally.

(2). Press YY for 1 second, the next program value to be set (current segment running time) will be displayed in turn, each segment of temperature control will be arranged in the order of SP T SP. That is the meaning as temperature of beginning of the segment--running time of the segment--target temperature value. The target value of temperature is the beginning temperature of next segments programs. Press XX after finished setting "-121.0" to quit the programs interface (Press XYZ to modify the value, after finished wanted value, press XX to confirm)

(3) Press YY for 2 seconds, user can return to the former parameter of programs.

(4) Press X firstly then press Y simultaneously to quit the temperature control setting state. Controller will quite the parameter setting state if without any operation for 10 seconds. (Please check and confirm the programs are set rightly after setting finished to start running of temperature controller.


temperature programs setting

2. Example of
temperature-time
programs setting

The controller adopts the form of SP


T SP to input the temperature curve,
and uses various symbol to prompt
the data to be input, and the shape of
the curve is determined by the
coordinates at the turning points.

Please fill out the data sheet in the following order and format before entering data:

Prompt	Input value	Meaning
SP 01	0	The beginning temperature
T- 01	20.0	Running time of first segment program
SP 02	200	Temperature value of the first turning point (Target temperature value of SP 01, The beginning temperature value of SP 02)
T- 02	120.0	Running time of SP 02
SP 03	200	Temperature value of the second turning point (Target temperature value of SP 02 , The beginning temperature value of SP 03)
T-03	80.0	Running time of SP 03
SP 04	1000	Temperature value of the third turning point (Target temperature value of SP 03 , The beginning temperature value of SP 04)
T-04	30.0	Running time of SP 04
SP 05	1000	Temperature value of the forth turning point (Target temperature value of SP 04 , The beginning temperature value of SP 05))
T- 05	100.0	Running time of SP 05
SP 06	500	Temperature value of the fifth turning point (Target temperature value of SP 05 , The beginning temperature value of SP 06)
T-06	-121.0 0	Finished programs setting


Using XXYY,input the temperature-time curve data into the controller in turn,then the programs setting finished.Attention:After finished the temperature curve,the end "-121.0" is a must to input!!!!When programs running,there is no access to modify the data.If in need to modify the programs,stop the programs is needed then modify the temperature curve.

4. Running of intelligent temperature controller

Automatic control is as below:

1. If the controller is only in the basis state(programs-stop,SV show STOP),Press XX for 1 second,enter running program state(PV-STEP,SV-XX segment),user can choose the needed step-segments by press ZZ or TT,Generally,the step number will raised automatically,and there is no need to interfere by operator.Sometimes,user wanted specified segments to start or jump at certain programs,then modify the STEP to achieve the goal.Then press XX+YY to return the basis state.
2. Press YY for 2seconds(SV show run),the controller enter into automatic control state.

5.HOLD(PAUSE) OF intelligent temperature controller

Pause is as below:

Press ZZ for 2 seconds when programs running,SV show “HOLD”symbol then the controller is at pause state(Temperature remain the control pause temperature and time stops increasing).In a paused state press TT again for 2 seconds ,SV show “RUN”,then controller will run again.

Program running state


Program running state Program continue to run”

Program
holding

6. Stop of Intelligent temperature controller

Stop control is as following:

When the programs is running or paused,press ZZ for 2 seconds,the controller will display “STOP”symbol,At this point,the program control is ended,the controller is in the basic state of stop state,meanwhile,STEP is modified to “1”.At this time,PV show furnace temperature”XXXX℃”, SV display“STOP”

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Maintenance precautions

Note:

After the power supply is turned on, the sound of the cooling fan running should be heard. If there is no such sound, stop the machine immediately to check and replace the cooling fan. Fan specification: Axial flow fan AC220v/120×120mm

If the furnace is put aside for a long time or works in a humid environment, the electric furnace should be baked at a low temperature for 3 hours when it is used for the first time [≤ 300 degrees]

Please place the fire block correctly when using the furnace.

1. When the furnace is in use, the input power must not exceed the rated power, and the furnace temperature must not exceed the rated temperature, so as not to damage the heating element and furnace lining.
2. The rising rate of each temperature section should not be too different. When setting the heating rate, the physical and chemical properties of the sintered material should be fully considered to avoid material spraying and pollute the furnace.
3. Regularly check whether the contact of the electrical connection part of the temperature control system is good.
4. This furnace is suitable for the following working conditions:
 - (6) The ambient temperature is between -10°C and 75°C .
 - (7) The relative humidity of the surrounding environment should not exceed 85%.
 - (8) There is no conductive dust, explosive gas and corrosive gas that seriously damages metal and insulating materials around the furnace.
 - (9) There is no obvious tilt, vibration and bump.
 - (10) When the electric furnace is put aside for a long time or the working environment is relatively humid, it is necessary to pay attention to the furnace drying operation for 3 hours [≤ 300 degrees] during use, so as to prevent the furnace from being damp and cracked, which will affect the heat preservation effect.